

# Polyethylene Glycol 1000


Application #AN3900

## Conditions

Column: ACE 3 C8  
Dimensions: 150 x 4.6 mm  
Part Number: ACE-112-1546  
Mobile Phase: A: 0.1% formic acid in H<sub>2</sub>O  
B: MeOH

Time (mins)	%B
0	50
45	85
50	50
60	50

Flow Rate: 1 mL/min  
Detection: APCI (negative ion)



Reproduced with permission of Wickham Laboratories, UK  
© 2015 ACT Ltd. All rights reserved. All products are available worldwide.