


# Determination of Baicalin in Children's Liquid Oral Medicine

**ACE**<sup>®</sup>  
Ultra-inert  
UHPLC & HPLC Columns

Application #AN6290

## Conditions

Column:	ACE Excel 5 AQ
Dimensions:	150 x 4.6 mm
Part Number:	EXL-126-1546U
Mobile Phase:	0.1% H <sub>3</sub> PO <sub>4</sub> in H <sub>2</sub> O/ MeOH (50:50 v/v)
Flow Rate:	1.0 mL/min
Injection:	5 µL
Temperature:	40 °C
Detection:	UV, 280 nm
Monograph Method:	ChP.2015, Vol. 1



1. Baicalin

