

Neomycin in HILIC

Conditions

Column: ACE 5 HILIC-A
Dimensions: 150 x 4.6 mm
Part Number: HILA-5-1546U
Mobile Phase: 20 mM ammonium formate (pH 4.7) in 80:20 MeCN/H₂O

Application #AN6930

Flow Rate: 1 mL/min
Injection Volume: 20 µL
Temperature: 30 °C
ELSD: Evaporator Temp: 50 °C
Gain: 4
Filter: 10 s
Sampling rate: 100 Hz

